

Basket Connection Guide

	SI Digital PC+	SISA / SISA 'S'	SI PCF / SI PCF 'S'	PCF PRO	Niagara
General Instruments	Basket & Hose Set: MED1018.1 Hose Set: MP22203	SISA Basket & Hose Set: MP21853 Hose Set: MP22203	Basket & Hose Set: MP21832 Hose Set: MP22183 Y Hose Set: MP22199	Basket & Hose Set: MP21852 Hose Set: MP21846	Carrier: MED1437 Baskets only: M20044 / M20043
		SISA 'S' Basket & Hose Set: MP21854 Hose Set: MP22203			Basket & Hose Set: MED1425 Hose Set: MP22243
					Basket: MP21841 Hose Set: MP22222
Orthopaedic Sets			Basket only: MED8603		Carrier: MED1415
SI EndoWrist® 8 mm			Basket & Hose Set: MP22062 Hose FP1: M20066	Basket & Hose Set: MP21907 Hose Set: MP22188	Carrier: MED1455 ⁹ Or Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21838 ¹ Hose Flush Port 1: MP22212
Si EndoWrist® 5 mm / SS			Basket & Hose Set: MP21833 Hose FP1: M20066 Hose for Shaft Port: MP22201	Basket & Hose Set: MP21907 Hose Set: MP22189	Carrier: MED1426 Hose Set: M20134
Si Stapler			Basket & Hose Set: MP21833 ⁷ Hose FP1: M20066 Hose for Shaft Port: MP18699.10	Basket & Hose Set: MP21907 Hose Set: MP22190	Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21839 ² Hose Flush Port 1: MP22214 Hose for Shaft Port: MP22215
Xi EndoWrist® 8 mm			Basket & Hose Set: MP21104 ⁹ Hose Set: MP21152	Basket & Hose Set: MP22036 Hose Set: MP22018	Carrier: MED1426 Hose Set: M20134 Hose Set: MP20966.5 Basket: MP21837 ³ Hose Set: MP21133 ⁴
Xi Stapler			Basket & Hose Set: MP21104 ⁸ Hose Set: MP21152	Basket & Hose Set: MP21906 Hose Set: MP22018	Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21840 ⁵ Hose Set: MP22217 Basket & Hose Set: MP21836 ⁶ Hose Set: MP21133 ⁴

1 Wire basket for 5 Si 8mm EndoWrists®

2 Sheet metal basket for 2 Si Staplers

3 Wire basket for 5 Xi EndoWrists®

4 Only to be used with MP21837 and MP21836 baskets

5 Sheet metal basket for 2 Xi Staplers

6 Wire basket for 3 Xi EndoWrists® / 2 Xi Staplers

7 Compatible with Si 8mm, Si 5mm, Si Single Site and Si Stapler

8 Basket including connector / hose for up to 4 Xi EndoWrists® and 2 Xi Staplers

9 Can only accommodate the Si 8mm EndoWrists® directly on carrier without the use of baskets