

Disposable Face Visors

The Valisafe Face Shields are designed to protect your eyes and face from accidents, by preventing exposure from hazardous or infectious substances. Perfect for all staff involved in equipment decontamination, medical, dental and emergency medical healthcare workers.

- Latex free
- Anti Fog (both sides)
- Protects from harmful debris & splashes
- Contoured foam for better fit and comfort
- Lightweight & comfortable to wear
- Fits comfortably with a surgical mask & over prescription eyewear



Standard length



Longer length



Standard Length:	Width: 33cm	Length: 19cm
Longer Length:	Width: 33cm	Length: 25cm

Order details:

V3505040 - Standard Length	Disposable Face Shields	200 faceshield frames
V3505040.1 - Longer Length	Disposable Face Shields	200 faceshield frames